

Ki67 Analyte Control^{DR}

Cell line controls for immunohistochemistry.

Research Use Only

PRODUCT AVAILABILITY

Product Code	Product Description
HCL065	X2 Cut slides
HCL066	X5 Cut slides
HCL067	X1 Cell microarray block

APPLICATION

This product is comprised of three human cell lines with low, intermediate and high expression of the nuclear proliferation marker Ki67. It is suitable for use as a same slide control material in immunohistochemistry. The control is intended for analytical assessment of reagent and instrument performance only.

The control can be used in combination with Visiopharm Qualitopix™ image analysis. Before use, every block should be validated to establish a two standard deviation QC threshold for each cell line in Qualitopix™. If any cell line falls out of the two standard deviation boundary, staining must be reviewed. If consecutive Qualitopix™ scores fall out of the boundary, the results of the test sample should be re-evaluated and the reagents and instrument investigated.


MATERIALS

Three formalin fixed paraffin embedded cell lines.


Cell line A: Urinary bladder carcinoma
Cell line B: Breast ductal carcinoma
Cell line C: Colorectal adenocarcinoma

Cells are fixed in 10% neutral buffered formalin and paraffin wax embedded.

Sections are cut at 4µm, mounted on positively charged slides and dried overnight at 37°C.



Cell microarrays (CMA) contain cores that are 1.5-2mm in diameter and 3-3.5mm in length. It is possible to obtain over 300 sections depending on thickness.



EXPRESSION PROFILE

Cell Line	IHC for Ki-67
A	Low
B	Intermediate
C	High

If any of the cell lines in the control fail to demonstrate the expected staining, results of the test sample should be re-evaluated.

STORAGE AND HANDLING

Store at 2-8°C. Do not freeze (for expiration date please see the product label)

WARNINGS AND PRECAUTIONS

1. The product is intended for research use only.
2. This product is not a diagnostic aid.
3. It is the responsibility of the end user to determine suitability with their reagents and procedures within their laboratory.
4. Do not use after expiration date printed on product labels. The user must validate any storage conditions other than those specified in the package insert.

TROUBLE SHOOTING

For further help please contact HistoCyte Laboratories Ltd at info@histocyte.com or call +44 (0)191 603 1007.



For updates and additional product information please visit: www.HistoCyte.com