

Product Data Sheet



Quality in Control

PRODUCT NAME: HPV/p16 Analyte Cells^{DR}


PRODUCT CODE: HCL001 (2 unstained slides)
HCL002 (5 unstained slides)

INTENDED USE: Research Use Only (RUO)

STORAGE: 2-8°C

DESCRIPTION: Each slide includes cell lines with a Dynamic Range (DR) of expression for a variety of biomarkers including HPV and p16. Each core has a 2mm diameter:

Cell line A: human breast adenocarcinoma
Cell line B: human cervical squamous cell carcinoma
Cell line C: human cervical adenocarcinoma
Cell line D: human epidermoid carcinoma



Fixative: 10% Neutral Buffered Formalin
Embedding: In paraffin wax
Section Thickness: 4µm
Mounting: Mounted on positively charged slides and dried at 37°C overnight