

Product Data Sheet



Quality in Control

PRODUCT NAME: HPV/p16 Analyte Cells


PRODUCT CODES: HCL004 (2 slide pack)
HCL005 (5 slide pack)

INTENDED USE: Research Use Only (RUO)

STORAGE: 2-8°C

DESCRIPTION: Each slide includes 3 cell lines of a 2mm diameter:

Cell line A: human breast adenocarcinoma
Cell line B: human cervical adenocarcinoma
Cell line C: human epidermoid carcinoma



Fixative: 10% Neutral Buffered Formalin
Embedding: In paraffin wax
Section Thickness: 4µm
Mounting: Mounted on positively charged slides and dried at 37°C overnight