

Product Data Sheet

PRODUCT NAME: HPV/p16 Analyte Cells

PRODUCT CODE: HCL006

INTENDED USE: Research Use Only (RUO).

STORAGE: 2°C to 8°C. Avoid freezing as this may cause the wax to crack.

DESCRIPTION: Cell microarray (CMA) wax block, which includes 3 individual cores (diameter 2.0mm x length 3.0-3.5mm) containing cell lines of the following tumour origin:

Cell line A: human breast adenocarcinoma
Cell line B: human cervical adenocarcinoma
Cell line C: human epidermoid carcinoma

Fixative: 10% Neutral Buffered Formalin
Embedding: In paraffin wax

DIRECTIONS FOR USE:

Material is supplied in a block (see figure 1 below) and sections should be cut and arranged as appropriate.




Figure 1. HPV/p16 Analyte Cells CMA block.