

Product Data Sheet



Quality in Control

PRODUCT NAME: ALK-Lung Analyte Cells

PRODUCT CODE: HCL007 (2 unstained slides)
HCL008 (5 unstained slides)


INTENDED USE: Research Use Only (RUO)

STORAGE: 2-8°C

DESCRIPTION: Each slide includes 2 cell lines of a 2mm diameter:

Cell line A: Breast adenocarcinoma

Cell line B: Pulmonary carcinoma with EML4-ALK gene translocation



Fixative: 10% Neutral Buffered Formalin

Embedding: In paraffin wax

Section Thickness: 4µm

Mounting: Mounted on positively charged slides and dried at 37°C overnight