

Product Data Sheet



Quality in Control

PRODUCT NAME: ALK-Lung Analyte Cell block

PRODUCT CODE: HCL009 (1 cell microarray (CMA) block)

INTENDED USE: Research Use Only (RUO)

STORAGE: 2°C to 8°C. Avoid freezing as this may cause the wax to crack

DESCRIPTION: CMA which includes 2 individual cores (Diameter 2.0mm x Length 3.0-3.5mm) containing cell lines of the following tumour origin:

Cell line A:	Breast adenocarcinoma
Cell line B:	Pulmonary carcinoma with EML4-ALK gene translocation

Fixative:	10% Neutral Buffered Formalin
Embedding:	In paraffin wax

DIRECTIONS FOR USE:

Material is supplied in a block (see figure 1 below) and sections should be cut and arranged as appropriate.




Figure 1. ALK-Lung Analyte Cells CMA block.