

Product Data Sheet

PRODUCT NAME: ALK-Lymphoma Analyte Cells


PRODUCT CODE: HCL010 (2 unstained slides)
HCL011 (5 unstained slides)

INTENDED USE: Research Use Only (RUO)

STORAGE: 2-8°C

DESCRIPTION: Each control slide includes 2 cell lines of a 2mm diameter:

Cell line A: Breast adenocarcinoma
Cell line B: Anaplastic large cell lymphoma with NPM-ALK gene translocation



Fixative: 10% Neutral Buffered Formalin
Embedding: In paraffin wax
Section Thickness: 4µm
Mounting: Mounted on positively charged slides and dried at 37°C overnight