

Product Data Sheet

PRODUCT NAME: Breast Analyte Cells


PRODUCT CODE: HCL013 (2 unstained slides)
HCL014 (5 unstained slides)

INTENDED USE: Research Use Only (RUO)

STORAGE: 2-8°C

DESCRIPTION: Each control slide includes 2 cell lines of a 2mm diameter:

Cell line A: Osteosarcoma
Cell line B: Breast ductal carcinoma



Fixative: 10% Neutral Buffered Formalin
Embedding: In paraffin wax
Section Thickness: 4µm
Mounting: Mounted on positively charged slides and dried at 37°C overnight