

# Product Data Sheet

**PRODUCT NAME:** Breast Analyte Cells

**PRODUCT CODE:** HCL015 (1 cell microarray (CMA) block)

**INTENDED USE:** Research Use Only (RUO)

**STORAGE:** 2°C to 8°C. Avoid freezing as this may cause the wax to crack

**DESCRIPTION:** CMA which includes 2 individual cores (Diameter 2.0mm x Length 3.0-3.5mm) containing cell lines of the following tumour origin:

Cell line A: Osteosarcoma  
Cell line B: Breast ductal carcinoma

Fixative: 10% Neutral Buffered Formalin  
Embedding: In paraffin wax

**DIRECTIONS FOR USE:**

Material is supplied in a block (see figure 1 below) and sections should be cut and arranged as appropriate.




Figure 1. Breast Analyte Cells CMA block.