

Product Data Sheet

PRODUCT NAME: Breast Analyte Cells^{DR}


PRODUCT CODE: HCL016 (2 unstained slides)
HCL017 (5 unstained slides)

INTENDED USE: Research Use Only (RUO)

STORAGE: 2-8°C

DESCRIPTION: Each slide includes 5 cell lines with a Dynamic Range (DR) of expression for a variety of biomarkers. The cell lines are of the following origin:

Cell line A:	Breast ductal carcinoma
Cell line B:	Breast ductal carcinoma
Cell line C:	Breast ductal carcinoma
Cell line D:	Breast adenocarcinoma
Cell line E:	Osteosarcoma



Fixative:	10% Neutral Buffered Formalin
Embedding:	In paraffin wax
Section Thickness:	4µm
Mounting:	Mounted on positively charged slides and dried at 37°C overnight