

Product Data Sheet



PRODUCT NAME: PD-L1 Analyte Cells^{DR}


PRODUCT CODE: HCL019 (2 unstained slides)
HCL020 (5 unstained slides)

INTENDED USE: Research Use Only (RUO)

STORAGE: 2-8°C

DESCRIPTION: Each slide includes 4 cell lines with a Dynamic Range (DR) of expression for PD-L1. Each core has a 2mm diameter:

Cell line A:	Breast ductal carcinoma
Cell line B:	Osteosarcoma
Cell line C:	Fibrosarcoma
Cell line D:	T cell non-Hodgkin Lymphoma



Fixative: 10% Neutral Buffered Formalin

Embedding: In paraffin wax

Section Thickness: 4µm

Mounting: Mounted on positively charged slides and dried at 37°C overnight