

# Product Data Sheet

**PRODUCT NAME:** ROS1 Analyte Cells

**PRODUCT CODE:** HCL022 (2 unstained slides)  
HCL023 (5 unstained slides)

**INTENDED USE:** Research Use Only (RUO)

**STORAGE:** 2-8°C

**DESCRIPTION:** Each control slide includes 2 control cell line cores of a 2mm diameter:

Cell line A: Breast adenocarcinoma

Cell line B: Pulmonary carcinoma with ROS1 gene translocation


Fixative: 10% Neutral Buffered Formalin

Embedding: In paraffin wax

Section Thickness: 4µm

Mounting: Mounted on positively charged slides and dried at 37°C overnight

**SLIDE LAYOUT:**



For more information, contact [info@histocyte.com](mailto:info@histocyte.com) or visit our website [www.histocyte.com](http://www.histocyte.com).