

# Product Data Sheet



Quality in Control

**PRODUCT NAME:** HER2 Analyte Cells<sup>DR</sup>


**PRODUCT CODE:** HCL026 (2 unstained cell microarray slides)  
HCL027 (5 unstained cell microarray slides)

**INTENDED USE:** Research Use Only (RUO)

**STORAGE:** 2°C to 8°C.

**DESCRIPTION:** Each slide has 4 individual cores (Diameter 2mm x Length 3-3.5mm) containing cell lines with a dynamic range of HER2:

Cell line A:	Breast adenocarcinoma
Cell line B:	Breast adenocarcinoma
Cell line C:	Gastric adenocarcinoma
Cell line D:	Breast adenocarcinoma



**Fixative:** 10% Neutral Buffered Formalin

**Embedding:** Paraffin wax

**Section thickness:** 4µm

**Mounting:** Mounted on positively charged slides and dried at 37°C overnight.