

# Product Data Sheet



Quality in Control

<b><u>PRODUCT NAME:</u></b>	HER2 Analyte Cells <sup>DR</sup>
<b><u>PRODUCT CODE:</u></b>	HCL028 (1 cell microarray (CMA) block)
<b><u>INTENDED USE:</u></b>	Research Use Only (RUO)
<b><u>STORAGE:</u></b>	2°C to 8°C. Avoid freezing as this may cause the wax to crack
<b><u>DESCRIPTION:</u></b>	CMA block that includes 4 individual cores (Diameter 2mm x Length 3-3.5mm) containing cell lines with a dynamic range of HER2:  Cell line A: Breast adenocarcinoma Cell line B: Breast adenocarcinoma Cell line C: Gastric adenocarcinoma Cell line D: Breast adenocarcinoma  Fixative: 10% Neutral Buffered Formalin Embedding: Paraffin wax

## **DIRECTIONS FOR USE:**

Material is supplied in a block (see figure 1) and should be cut and arranged as appropriate.




Figure 1. HER2 Analyte Cell CMA block.