

HCL053_HCL054_HCL055_PI_V001

ALK Analyte Control^{DR}

Cell line controls for immunohistochemistry and in situ hybridization.

Research Use Only

PRODUCT AVAILABILITY

Product Code	Product Description	
HCL053	X2 Cut slides	
HCL054	X5 Cut slides	
HCL055	X1 Cell microarray block	

APPLICATION

This product is suitable for use in immunohistochemistry and in situ hybridization.

MATERIALS

Four formalin fixed paraffin embedded cell lines.

Cell line A: Breast adenocarcinoma Cell line B: Glioblastoma Cell line C: Lung adenocarcinoma Cell line D: Anaplastic large cell lymphoma

Cells are fixed in 10% neutral buffered formalin and paraffin wax embedded.

Sections are cut at $4\mu m,$ mounted on positively charged slides and dried overnight at 37°C.



Cell microarrays (CMA) contain cores that are 1.5-2mm in diameter and 3-3.5mm in length. It is possible to obtain over 300 sections depending on thickness.



Expression Profile

Cell Line	IHC for ALK-EML4 (Lung)	IHC for ALK- NPM (Lymphoma)	FISH for ALK
Α	Negative	Negative	Negative
В	Low	Negative	Negative
С	High	Mid	Positive
D	High	High	Positive

Storage and Handling

Store at 2-8°C. Do not freeze (for expiration date please see the product label)

WARNINGS AND PRECAUTIONS

- 1. The product is intended for research use only.
- 2. It is the responsibility of the end user to determine suitability with their reagents and procedures within their laboratory.
- 3. Do not use after expiration date printed on product labels. The user must validate any storage conditions other than those specified in the package insert.

TROUBLE SHOOTING

For further help please contact HistoCyte Laboratories Ltd at info@histocyte.com or call +44 (0)191 603 1007.



For updates and additional product information please visit: www.HistoCyte.com